MARROW-PACK

EN	Instructions	P. 02
IT	Istruzioni	P. 06
EL	Οδηγίες	P. 10
PT	Instruções	P. 14
ES	Instrucciones	P. 18

Intended Use: the needle MARROW-PACK is a Bone marrow harvest needle designed for autologous bone marrow harvest transplantation procedures.


CHARACTERISTICS AND PERFORMANCE:

MARROW-PACK is a single-use device designed and manufactured to get a large volume of bone marrow. The needle is composed of a cannula and an internal stylet in stainless steel and a plastic handle with a female standard Luer-lock connection. The holes on the distal portion of the cannula are placed offset each other, in order to improve the explant of bone marrow.

MARROW-PACK features include:

- Ring nut and spacer that allows millimetric rising of the needle
- indication of the cm on the graduated scale
- holes on the distal cannula wall
- closed distal trocar tip

Kit includes an aspiration syringe with a standard male Luer-lock connection, a filtering system (to retain any aspirated bone fragments) and an injection syringe with a standard male Luer-lock connection.

This device has not to be used in combination with other medical devices.

INDICATIONS:

The device is indicated for bone marrow explant.

CONTRA-INDICATIONS:

Contraindicated for use on patients who are receiving heavy anticoagulant therapy or who have a severe bleeding problem. Physician judgement is required when considering biopsy on patients with a bleeding disorder or receiving anticoagulant medications.

Target Patient:

The device is intended for the following patient population:

SEX: Male – Female

EXCLUSION: None (the choice of diameter and length is up to the surgeon based on the age, feature, and pathologic conditions of the patient)

The device is used to treat patients with suspicion of disease (e.g. inflammation or tumour), to arrive at a histopathological diagnosis, and therefore to precisely define its characteristics (severity, extent, possible therapies).

The choice of the patient and the decision to perform a biopsy must be made by medical staff taking into account the full clinical and physical picture of the patient. The choice must be made by evaluating various aspects, for example, weight/fat mass of the patient, type of patient and any other physical/clinical characteristic of the patient that may have an impact on the biopsy result. Based on these data the expert and the trained surgeon may decide if perform the biopsy on the patient and which needle is more suitable to do it.

INTENDED USERS AND ENVIRONMENT USE CONDITIONS:

The device is for professional use and, in particular, must be used by qualified and trained surgeons.

This product should be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of biopsy needle use, in particular, those relating to the specific organ being biopsied.

The procedure has to comply with the antisepsis rules.

WARNINGS, PRECAUTIONS AND SIDE EFFECTS:


- Do not use the instrument if the package is damaged.
- Do not use the instrument beyond the expiry date.
- The device is single-use and must be destroyed after use. Reuse is strictly prohibited. Reuse or resterilization may compromise the structural integrity of the device and/or cause damage to the device which may cause injury, disease or death to the patient. Reuse or resterilization may also create a risk of device contamination and/or cause infection of the patient or cross-infection, including but not limited to transmission of infectious diseases from one patient to the next. Device contamination may cause harm, disease or death to the patient.
- It is also recommended not to use the device if particularly risky manoeuvres are necessary where it will be subjected to high torsional or bending forces which might compromise the integrity and proper functioning of the system.
- Post-biopsy patient care may vary with the biopsy technique utilized and the individual patient's physiological conditions. Observation of vital signs and other usual precautions to avoid the serious potential complications associated with any biopsy procedure may be indicated.
- Good medical judgement should be exercised in considering biopsy on patients who are receiving anticoagulant therapy or who have bleeding disorders.
- This product should be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of needle biopsy, in particular, those relating to the specific organ being biopsied.
- The choice of the needle type, gauge and size to be used must be made by medical staff taking into account the full clinical and physical picture of the patient. The choice must be made by evaluating various aspects. For example weight/fat mass of the patient, type of patient and any other physical/clinical characteristic of the patient that may have an impact on the needle type to be used.

INSTRUCTION FOR USE**Preliminary operations**


- Before use check that the package is not damaged and that the set and its components are not damaged or kinked; if so, inform the manufacturer or your area representative. Take the devices out of the package operating in aseptic conditions.
- Sedate and position the patient.
- Spray the instruments to be used during the suction procedure with anticoagulant and fill the suction syringe with an amount of anticoagulant at a ratio of 1:5 with respect to the volume of tissue to be suctioned.
- Perform the necessary skin disinfection procedures.

Instructions:


After positioning and preparing the patient, insert the needle perpendicular to the bone surface holding the grip resting on the palm of your hand.


Move the needle forward by turning it to the right and left at the same time pushing until you find the least resistance (which indicates that the needle has entered the marrow cavity).

Extract the internal stylet and draw the marrow tissue by inserting the dedicated aspiration syringe (the one with the rotating knob on the rear) into the needle Luer-lock cone and proceed with the first aspiration.


Reinsert the stylet into the needle. Then continue the forward movement down to the desired depth. Turn the rotating ring nut until the base of the depth limiter rests on the patient's skin. Proceed with the aspiration.



Perform millimetric needle retraction by turning the rotating ring nut and always holding the needle depth limiter resting on the patient's skin. Once the desired retraction point has been reached, proceed with a new aspiration.

Repeat the retraction and aspiration operation at the various levels desired.

Finally, remove the needle from the patient.

Detach the aspiration syringe and connect it to the filter in correspondence to the inlet connection. Connect the filter outlet connection to the second syringe provided in the kit (which must be fully closed). Taking care to keep the system so connected vertically with the suction syringe upwards, proceed with the infusion of the tissue aspirated by the first syringe and its consequent aspiration into the second syringe to filter the tissue.

PACKAGING, STERILIZATION, SINGLE-USE:

The system is individually packed in rigid PVC blister and medical-grade paper. Do not use the device if the package is open or damaged or contaminated. It is sold in packs of 1 piece.

The needle is sterilized with Ethylene Oxide (EtO) in compliance with ISO 11135 requirements.

The device is not reusable. The device is single-use and must be destroyed after use.

DISPOSAL:

After use, the device may constitute a risk of contamination and/or infection and puncture or injuries caused by sharp parts of the device. Therefore, handle with care and dispose of the product in accordance with the regulations in force and medical practice.

NOTICE FOR USER AND/OR PATIENT: ANY SERIOUS INCIDENT THAT HAS OCCURRED IN RELATION TO THE DEVICE SHOULD BE REPORTED TO THE MANUFACTURER AND THE COMPETENT AUTHORITY OF THE MEMBER STATE IN WHICH THE USER AND/OR PATIENT IS ESTABLISHED.

Biopsybell declines all responsibility for any damage caused by improper use or different from what is indicated in this instruction leaflet.

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